## DO-IT-YOURSELF BACTERIA IN WATER TEST

Our (DIY) Safe Home® Bacteria in Water Test vials are designed to be used in remote areas or in the comfort of your own home. Simply follow the detailed instructions on the reverse side of this paper.

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## CONTENTS NEEDED FOR THIS TEST

One (DIY) test vial with white "tamper evident" cap (contains test media)  
**These Instructions**

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## IMPORTANT INFORMATION

**DIY test media** - Composition: a mixture of organic and inorganic compounds. Ingredients include: Sodium Chloride, Tryptose, Dipotassium Hydrogenorthophosphate, Potassium dihydrogenorthophosphate, L-Tryptophan, D-Glucitol. Does not contain any chemicals known to cause cancer, birth, or any other reproductive defects.

**Handling DIY Test Vial** - Before use, contents of this vial should be stored, tightly closed, dry, protected from light and at a temperature between 59° F to 77° F (15° C to 25° C). Keep away from strong oxidizing agents. Do not touch the mouth of the vial or the inside of the cap. You will view the test results in the comfort of your own home. Do not return this DIY test vial to our laboratory. Once test is completed, dispose of test vial and its contents as directed. Keep out of reach of children.

**First Aid** - Do not ingest test media. If test media is inhaled, seek fresh air. If test media makes skin contact, wash off with plenty of water and remove contaminated clothing. If test media makes eye contact, rinse eyes out with plenty of water. If test media is swallowed, drink water (two glasses at most). Consult a doctor if feeling unwell.

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## INSTRUCTIONS FOR DO-IT-YOURSELF BACTERIA IN DRINKING WATER SAMPLING AND TESTING

## SAMPLE COLLECTION

1. Wash your hands before starting. Remove the test vial with the white "tamper evident" cap.
2. Turn on cold water and allow it to run for approximately five minutes before sampling. Reduce the flow of water to about the diameter of a pencil. Remove the white cap completely. Tilt the vial slightly and slowly fill the vial to the 25 mL mark near the top. **Do not fill past this mark. Do not overflow the vial.**

3. Tightly screw the cap on the vial and shake well, making sure that the media in the vial dissolves completely. Set the vial in an upright position and on a flat surface.

## SAMPLE INCUBATION

4. Incubate your water sample at temperatures in the following ranges (70°-105° F or 21°-40.6° C) indoors or outdoors. A temperature of 95.5° F / 35.3° C is best.


Test results as fast as \*24 hours (\*speed of test results may vary). Test results for bacteria will take a minimum of 24 hours (± 4 hours) to complete. Lower incubation temperatures and/or lower bacteria numbers may cause final results to take up to 72 hours.

Store vial in an upright position (out of direct sunlight) until test's result is evident.

Shaking the vial for 5-10 seconds (every 8-12 hours) can enhance color development but is not required.

## TEST RESULTS

5. After proper incubation time has elapsed and results are evident, without opening the vial, observe the color of the liquid.

**Yellow = Absent** (negative result for coliform bacteria). Coliform bacteria, including E. coli bacteria, are absent. Note: As long as the sample inside the vial is clear or yellow, allow incubation to continue for at least 72 hours to confirm the absence of bacteria.

**Teal or Blue/Green = Present** (positive result for coliform bacteria). This test result could include E. coli. Color and intensity may vary. **NOTE: If any shade of blue/green exists at any time during the test period (between 24-72 hours), the test is complete.**

6. When the test is complete, open vial and discard liquid into toilet. Place a few drops of household bleach into the vial, cap it, shake it, and once again discard the contents into the toilet. Recap the vial. Place the vial in a trash can. Wash your hands thoroughly with soap and warm water.

The results of your **Safe Home (DIY) Bacteria in Water Test** are quite accurate. This test is a "Presence/Absence" indicator of bacteria. If you desire additional testing, simply go to [safehometestkits.com](http://safehometestkits.com) and choose from the variety of drinking water test kits offered on our website.

(DIY Bacteria Test Instructions for Single Test and Multi-Packs)

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